



# Proteogenix

**Proteogenix SAS**

www.proteogenix.science  
19 rue de la Haye  
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70  
Fax +33 (0) 978 53 36 90  
products@proteogenix.fr

## Anti-MERS-CoV RBD Monoclonal Antibody, AT2F7

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PTX19385
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG2a, kappa
<b>TARGET</b>	Spike glycoprotein, S glycoprotein, E2, Peplomer protein, RBD, Spike protein S1, S, Receptor-Binding Domain, SARS-like Coronaviruses

### PRODUCT INFORMATION

---

<b>HOST</b>	Mouse
<b>REACTIVITY</b>	Middle East respiratory syndrome-related coronavirus (isolate United Kingdom/H123990006/2012) (MERS-CoV) (Betacoronavirus England 1)
<b>FORM &amp; BUFFER</b>	0.01M PBS, pH 7.4.
<b>STORAGE</b>	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



**WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.**