



Anti-MERS-CoV RBD Monoclonal Antibody, JC57-14

OVERVIEW

PRODUCT REFERENCE	PTX19386
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	Spike glycoprotein, S glycoprotein, E2, Peplomer protein, RBD, Spike protein S1, S, Receptor-Binding Domain, SARS-like Coronaviruses

PRODUCT INFORMATION

HOST	Human
REACTIVITY	Middle East respiratory syndrome-related coronavirus (isolate UnitedKingdom/H123990006/2012) (MERS-CoV) (Betacoronavirus England 1)
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

