



## Anti-MERS-CoV RBD Monoclonal Antibody, CDC2-C2

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PTX19387
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1
<b>TARGET</b>	Spike glycoprotein, S glycoprotein, E2, Peplomer protein, RBD, Spike protein S1, S, Receptor-Binding Domain, SARS-like Coronaviruses

### PRODUCT INFORMATION

---

<b>HOST</b>	Human
<b>REACTIVITY</b>	Middle East respiratory syndrome-related coronavirus (isolate UnitedKingdom/H123990006/2012) (MERS-CoV) (Betacoronavirus England 1)
<b>FORM &amp; BUFFER</b>	0.01M PBS, pH 7.4.
<b>STORAGE</b>	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

