



Anti-MERS-CoV S2 Protein Monoclonal Antibody, G4

OVERVIEW

PRODUCT REFERENCE	PTX19389
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	Spike glycoprotein, S glycoprotein, E2, Peplomer protein, Spike protein S2, S, S2 Protein, SARS-like Coronaviruses

PRODUCT INFORMATION

HOST	Human
REACTIVITY	Middle East respiratory syndrome-related coronavirus (isolate UnitedKingdom/H123990006/2012) (MERS-CoV) (Betacoronavirus England 1)
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

