



Anti-MARV GP/Envelope glycoprotein Monoclonal Antibody, MR191

OVERVIEW

PRODUCT REFERENCE	PTX19398
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	Envelope glycoprotein, GP1, 2, GP, Virion spike glycoprotein, GP1, GP2, GP, Marburg Virus (MARV)

PRODUCT INFORMATION

HOST	Human
REACTIVITY	Lake Victoria marburgvirus (strain Ravn-87) (MARV) (Marburg virus (strainKenya/Ravn/1987))
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

