



# Proteogenix

**Proteogenix SAS**

www.proteogenix.science  
19 rue de la Haye  
67300 Schiltigheim, France  
Tel +33 (0) 390 20 54 70  
Fax +33 (0) 978 53 36 90  
products@proteogenix.fr

## Anti-MARV GP/Envelope glycoprotein Monoclonal Antibody, MR78

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PTX19399
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1
<b>TARGET</b>	Envelope glycoprotein, GP1, 2, GP, Virion spike glycoprotein, GP1, GP2, GP, Marburg Virus (MARV)

### PRODUCT INFORMATION

---

<b>HOST</b>	Human
<b>REACTIVITY</b>	Lake Victoria marburgvirus (strain Ravn-87) (MARV) (Marburg virus (strainKenya/Ravn/1987))
<b>FORM &amp; BUFFER</b>	0.01M PBS, pH 7.4.
<b>STORAGE</b>	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



**WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.**