



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-MARV VP35/Marburg VP35 Monoclonal Antibody, SAA1444

OVERVIEW

PRODUCT REFERENCE	PTX19401
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	Polymerase cofactor VP35, Marburg VP35, mVP35, VP35, Marburg Virus (MARV)

PRODUCT INFORMATION

HOST	Human
REACTIVITY	Lake Victoria marburgvirus (strain Popp-67) (MARV) (Marburg virus (strain West Germany/Popp/1967))
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.