



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-Human CD95/FAS Monoclonal Antibody, F919

OVERVIEW

PRODUCT REFERENCE	PTX19402
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	APT1, TNFRSF6, FAS1, FASLG receptor, CD95, Apo-1 antigen, Tumor necrosis factor receptor superfamily member 6, FAS, Apoptosis-mediating surface antigen FAS

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.