



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-Human MYH2 Monoclonal Antibody, SF9

OVERVIEW

PRODUCT REFERENCE	PTX19403
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a
TARGET	Myosin heavy chain 2, MYHSA2, Myosin heavy chain IIa, MyHC-2a, Myosin-2, Myosin heavy chain, skeletal muscle, adult 2, MYH2, MyHC-IIa, Myosin heavy chain 2a

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.