



Anti-Human MYH9 Monoclonal Antibody, SAA1446

OVERVIEW

PRODUCT REFERENCE	PTX19406
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a
TARGET	NMMHC-IIA, Non-muscle myosin heavy chain IIa, Cellular myosin heavy chain, type A, Myosin heavy chain 9, NMMHC II-a, Non-muscle myosin heavy chain A, Myosin-9, MYH9, NMMHC-A, Myosin heavy chain, non-muscle IIa

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

