



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-Human GIPR/GIP-R Monoclonal Antibody, SAA1449

OVERVIEW

| | |
|--------------------------|---|
| PRODUCT REFERENCE | PTX19409 |
| CLONALITY | Monoclonal Antibody |
| ISOTYPE | IgG2a |
| TARGET | Gastric inhibitory polypeptide receptor, GIPR, GIP-R, Glucose-dependent insulinotropic polypeptide receptor |

PRODUCT INFORMATION

| | |
|--------------------------|---|
| HOST | Mouse |
| REACTIVITY | Human |
| FORM & BUFFER | 0.01M PBS, pH 7.4. |
| STORAGE | Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles. |



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.