



Anti-pan-MAGUK Monoclonal Antibody, K28/86

OVERVIEW

PRODUCT REFERENCE	PTX19410
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	PSD-95, PSD95, Disks large homolog 4, Synapse-associated protein 90, DLG4, Postsynaptic density protein 95, SAP-90, SAP90

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human, Mouse, Rat, Drosophila
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

