



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

VHH-hFc (IgG1) Isotype Control antibody (cAb-Lys2)

OVERVIEW

PRODUCT REFERENCE	PTX19420
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-hFc

PRODUCT INFORMATION

EXPRESSION SYSTEM

BUFFER 0.01M PBS pH 7.4 and 0.09% sodium azide

STORAGE CONDITIONS Store at 4 °C for short term. Freeze and store at -80 °C for long term.
Addition of 20 to 40% glycerol may improve cryoprotection.
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.
Aliquot the product according to your needs to avoid freeze/thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.