



## Anti-Human IL20RB Monoclonal Antibody

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PTX19430
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG2a, kappa
<b>TARGET</b>	IL-20R-beta, IL-20R2, FNDC6, IL-20RB, DIRS1, IL20RB, Interleukin-20 receptor subunit beta, Fibronectin type III domain containing 6, IL-20 receptor subunit beta

### PRODUCT INFORMATION

---

<b>HOST</b>	Mouse
<b>REACTIVITY</b>	Human
<b>FORM &amp; BUFFER</b>	0.01M PBS, pH 7.4.
<b>STORAGE</b>	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

