



Anti-Human IL20RB Polyclonal Antibody

OVERVIEW

PRODUCT REFERENCE	PTX19431
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human IL20RB (Ala37-Glu290)(PX-P6279)

PRODUCT INFORMATION

HOST	Rabbit
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4, 50% glycerol, 0.05% Proclin 300.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

