



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-Human IL24 Polyclonal Antibody

OVERVIEW

PRODUCT REFERENCE	PTX19432
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	E. coli - derived recombinant Human IL24 (Gln52-Leu206)(PX-P6280)

PRODUCT INFORMATION

HOST	Rabbit
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4, 50% glycerol, 0.05% Proclin 300.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.