



# Proteogenix

**Proteogenix SAS**

www.proteogenix.science  
19 rue de la Haye  
67300 Schiltigheim, France  
Tel +33 (0) 390 20 54 70  
Fax +33 (0) 978 53 36 90  
products@proteogenix.fr

## Anti-CD25/IL2RA Polyclonal Antibody

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PTX20637
<b>CLONALITY</b>	Polyclonal Antibody
<b>ISOTYPE</b>	IgG

### PRODUCT INFORMATION

---

<b>PRODUCTION HOST</b>	Rabbit
<b>SPECIES REACTIVITY</b>	Human
<b>FORM</b>	Liquid
<b>BUFFER</b>	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300
<b>TARGET</b>	Interleukin-2 receptor subunit alpha, IL-2 receptor subunit alpha, IL-2-RA, IL-2R subunit alpha
<b>STABILITY &amp; STORAGE</b>	Store at 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.



**WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.**