



Human Ade2 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1066
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Full-length
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	48,40 kDa
FORM	liquid
BUFFER	PBS, 20% glycerol, 0,05% sodium azide
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P22234
ALIAS/SYNONYMS	Ade2, multifunctional protein ADE2 isoform 2

