



Leishmania CPB Recombinant Protein

OVERVIEW

| | |
|---------------------------|------------|
| PRODUCT REFERENCE | PX-P1082 |
| ORIGIN SPECIES | Leishmania |
| FRAGMENT | Partial |
| TAG | Yes |
| DELIVERY CONDITION | Dry Ice |

PRODUCT INFORMATION

| | |
|--------------------------------|--|
| EXPRESSION SYSTEM | Prokaryotic expression |
| HOST | Escherichia coli (E.coli) |
| MOLECULAR WEIGHT | 18,82 kDa |
| FORM | liquid |
| BUFFER | PBS, Urea 8M |
| STABILITY & STORAGE | 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection) |
| UNIPROT | A4HTP0 |
| ALIAS/SYNONYMS | CPB, Cathepsin L-like Protease |

