



Trichinella spiralis GST24 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1104
ORIGIN SPECIES	Trichinella spiralis
FRAGMENT	Full-length
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	25,32 kDa
FORM	liquid
BUFFER	PBS, imidazole 300mM
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q0H8U8
ALIAS/SYNONYMS	GST24, glutathione-S-transferase

