



Mycoplasma LppA Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1122
ORIGIN SPECIES	Mycoplasma
FRAGMENT	Full-length
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	62,97 kDa
FORM	liquid
BUFFER	PBS, imidazole 250mM
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P55801
ALIAS/SYNONYMS	LppA, immunodominant protein P72, MSC_0014, P72

