



Trichinella spiralis NBL2 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1129
ORIGIN SPECIES	Trichinella spiralis
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	45,20 kDa
FORM	liquid
BUFFER	NaH ₂ PO ₄ 100mM, TrisHCl 10mM, Urea 8M
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	B0F9T2
ALIAS/SYNONYMS	NBL2ΔNter, hypothetical protein

