



Plasmodium falciparum RESA-2 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1146
ORIGIN SPECIES	Plasmodium falciparum
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	47,06 kDa
FORM	liquid
BUFFER	TrisHC 50mMl, 2% SDS, pH8
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q8IHM8
ALIAS/SYNONYMS	RESA-2, RESA-like protein with PHIST and DnaJ domains

