



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

S. pombe Smc1Nter Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1152
ORIGIN SPECIES	S. pombe
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	73,78 kDa
FORM	liquid
BUFFER	NaH ₂ PO ₄ 100mM, TrisHCl 10mM, Urea 8M
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	O94383
ALIAS/SYNONYMS	Smc1Nter, Structural maintenance of chromosomes protein 1



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.