



Staphylococcus AmiCat Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P2062
ORIGIN SPECIES	Staphylococcus
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	28.11 kDa
FORM	Frozen
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	KFA44079.1
ALIAS/SYNONYMS	AmiCat, chromogranin A (parathyroid secretory protein 1)

