



P. americana BetaGluc Recombinant Protein (E.coli production)

OVERVIEW

PRODUCT REFERENCE	PX-P2065
ORIGIN SPECIES	P. americana
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	58.76 kDa
FORM	Lyophilized
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	A0A059WHP2
ALIAS/SYNONYMS	BetaGluc, Beta-glucosidase

