



Dermatophagoides farinae (American house dust mite) Der f 1 allergen(Derf1) Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P3018
ORIGIN SPECIES	Dermatophagoides farinae (American house dust mite)
FRAGMENT	Full-length
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Insect
MOLECULAR WEIGHT	51.58 kDa
FORM	Frozen
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	A1YW12
ALIAS/SYNONYMS	Der f 1 allergen(Derf1)

