



Human CD81 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P3039
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	10.59 kDa
FORM	Frozen
BUFFER	PBS,pH 7.5 urea +8M
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P60033
ALIAS/SYNONYMS	CD81, S5.7, TAPA1, TSPAN28, Tetraspanin 28, Target Of Antiproliferative Antibody 1, 26 kDa cell surface protein TAPA-1

