



## Listeria monocytogenes Invasion associated protein p60 Recombinant Protein

### OVERVIEW

---

|                           |                        |
|---------------------------|------------------------|
| <b>PRODUCT REFERENCE</b>  | PX-P3056               |
| <b>ORIGIN SPECIES</b>     | Listeria monocytogenes |
| <b>FRAGMENT</b>           | Partial                |
| <b>TAG</b>                | Yes                    |
| <b>DELIVERY CONDITION</b> | Dry Ice                |

### PRODUCT INFORMATION

---

|                                |  |
|--------------------------------|--|
| <b>EXPRESSION SYSTEM</b>       | Prokaryotic expression   |
| <b>HOST</b>                    | Escherichia coli (E.coli)  |
| <b>MOLECULAR WEIGHT</b>        | 48.7 kDa   |
| <b>FORM</b>                    | Frozen   |
| <b>BUFFER</b>                  | PBS,pH 7.5   |
| <b>STABILITY &amp; STORAGE</b> | 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection) |
| <b>UNIPROT</b>                 | Q84DP3   |
| <b>ALIAS/SYNONYMS</b>          | iap, Invasion associated protein p60   |

