



Human Engrailed1 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P3062
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Full-length
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	42.14 kDa
FORM	Frozen
BUFFER	PO ₄ (pH:7.4) 20 mM with about 0.5M NaCl
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	
ALIAS/SYNONYMS	Engrailed1

