



Rainbow trout rtCD64s Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P3064
ORIGIN SPECIES	Rainbow trout
FRAGMENT	Unknown
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	34.74 kDa
FORM	Frozen
BUFFER	PBS pH 7.5 +urea 8M
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	
ALIAS/SYNONYMS	rtCD64s

