



Moloney murine leukemia virus M-MLV Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P3069
ORIGIN SPECIES	Moloney murine leukemia virus
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	70.26 kDa
FORM	Frozen
BUFFER	Tris-HCl 20mM(pH 7.5 at 25°C), 0.2M NaCl, 0.EDTA 1mM, DTT 1mM, 0.01% Nonidet P-40 and 50% glycerol.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P03355
ALIAS/SYNONYMS	M-MLV

