



## Human T-cell surface glycoprotein CD3 epsilon chain (CD3E) recombinant protein

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P4014
<b>ORIGIN SPECIES</b>	Homo sapiens (Human)
<b>FRAGMENT</b>	Partial
<b>TAG</b>	Yes
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Eukaryotic expression
<b>HOST</b>	Mammalian cells
<b>MOLECULAR WEIGHT</b>	15.09kDa
<b>FORM</b>	Lyophilized
<b>BUFFER</b>	50 mM Tris-HCl, pH 8.0, 150 mM NaCl
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	P07766
<b>ALIAS/SYNONYMS</b>	CD3E, T3E, TCRE, CD3-E Molecule, Epsilon(CD3-TCR Complex), T-cell surface antigen T3/Leu-4 epsilon chain

