



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Human Myoglobin recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P4024
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Full-length
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	18.06kDa
FORM	Lyophilized
BUFFER	50 mM Tris-HCl pH 7.5, 150 mM NaCl
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P02144
ALIAS/SYNONYMS	myoglobin, PVALB, MB, MGC13548



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.