



Human Beta-nerve Growth Factor proteins recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P4027
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Full-length
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	53.31kDa
FORM	Lyophilized
BUFFER	0.1M Tris-HCl, 0.1M Glycine, pH=7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P01138
ALIAS/SYNONYMS	NGF, NGFB, Beta-NGF, HSN5, NGF-B, Beta-Nerve Growth Factor proteins

