



Human AMH, Muellerian-inhibiting factor recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P4055
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Full-length
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	85.45kDa
FORM	Lyophilized
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P03971
ALIAS/SYNONYMS	MIF, MIH, MIS, Müllerian Inhibiting Factor, Müllerian Inhibiting Hormone, Müllerian Inhibiting Substance, AMH

