



## Rhesus monkey IL37 - Interleukin 1 Zeta (IL1z) recombinant protein

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P4063
<b>ORIGIN SPECIES</b>	Rhesus monkey
<b>FRAGMENT</b>	Full-length
<b>TAG</b>	Yes
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Prokaryotic expression
<b>HOST</b>	Escherichia coli (E.coli)
<b>MOLECULAR WEIGHT</b>	27.71kDa
<b>FORM</b>	
<b>BUFFER</b>	
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	
<b>ALIAS/SYNONYMS</b>	IL1F7, IL37, FIL1, FIL1(ZETA), FIL1Z, IL-1F7, IL-1H4, IL-1RP1, IL1H4, IL1RP1, Interleukin 37, Interleukin-1-Related Protein, Interleukin 1 Family, Member 7

