



## CD9 Recombinant Protein

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P4076
<b>ORIGIN SPECIES</b>	Homo sapiens (Human)
<b>FRAGMENT</b>	Ser112~Ile195
<b>TAG</b>	C-terminal His Tag
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Prokaryotic expression
<b>HOST</b>	Escherichia coli (E.coli)
<b>MOLECULAR WEIGHT</b>	10.5kDa
<b>FORM</b>	Frozen
<b>BUFFER</b>	PBS,pH 7.5
<b>STABILITY &amp; STORAGE</b>	Store at - 20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage.
<b>UNIPROT</b>	P21926
<b>ALIAS/SYNONYMS</b>	MIC3, TSPAN29

