



CD54 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P4085
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1~Glu480
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	53kDa
FORM	Frozen
BUFFER	PBS, pH7.5
STABILITY & STORAGE	Store at - 20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage.
UNIPROT	P05362
ALIAS/SYNONYMS	ICAM-1

