



CD98 Recombinant Protein

OVERVIEW

| | |
|---------------------------|----------------------|
| PRODUCT REFERENCE | PX-P4096 |
| ORIGIN SPECIES | Homo sapiens (Human) |
| FRAGMENT | Arg206-Ala630 |
| TAG | C-terminal His Tag |
| DELIVERY CONDITION | Dry Ice |

PRODUCT INFORMATION

| | |
|--------------------------------|---|
| EXPRESSION SYSTEM | Eukaryotic expression |
| HOST | Mammalian cells |
| MOLECULAR WEIGHT | 48kDa |
| FORM | |
| BUFFER | PBS pH7.5 |
| STABILITY & STORAGE | Store at - 20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. |
| UNIPROT | P08195 |
| ALIAS/SYNONYMS | SLC3A2,MDU1 |

