



CD121b Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P4102
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Glu343
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	45-50 kDa
FORM	Frozen
BUFFER	PBS, pH7.5
STABILITY & STORAGE	Store at - 20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage.
UNIPROT	P27930
ALIAS/SYNONYMS	IL1R2,IL-1R-2, IL-1RT2, CDw121b, IL-1R-beta,

