



CD257 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P4112
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Ala134–Leu285
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	17.97kDa
FORM	
BUFFER	PBS pH7.5
STABILITY & STORAGE	Store at - 20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage.
UNIPROT	Q9Y275
ALIAS/SYNONYMS	BAFF, BLYS, TALL1, TNFSF20, ZTNF4

