



## CD344 Recombinant Protein

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P4122
<b>ORIGIN SPECIES</b>	Homo sapiens (Human)
<b>FRAGMENT</b>	Met1-Glu180
<b>TAG</b>	C-terminal His Tag
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Eukaryotic expression
<b>HOST</b>	Mammalian cells
<b>MOLECULAR WEIGHT</b>	22-28kDa
<b>FORM</b>	
<b>BUFFER</b>	PBS pH7.5
<b>STABILITY &amp; STORAGE</b>	Store at - 20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage.
<b>UNIPROT</b>	Q9ULV1
<b>ALIAS/SYNONYMS</b>	EVR1,FEVR,Fz-4,Fz4,FZD4S,FzE4,GPCR,hFz4

