



Placenta-specific protein 1(PLAC1)

OVERVIEW

PRODUCT REFERENCE	PX-P4344
ORIGIN SPECIES	
FRAGMENT	Gln23-Met212
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	48.75kDa
FORM	
BUFFER	PBS pH7.5, 4M urea
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9HBJ0
ALIAS/SYNONYMS	

