



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Outer membrane virulence protein YopE (yopE)

OVERVIEW

PRODUCT REFERENCE	PX-P4378
ORIGIN SPECIES	
FRAGMENT	Met12-Met230
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	24.14kDa
FORM	
BUFFER	20 mM Tris-HCl pH8.0, 100 mM NaCl, 1 mM DTT, 0,1% sodium azide, 0.1mM EDTA, 10% glycerol
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P31492
ALIAS/SYNONYMS	Outer membrane virulence protein YopE, YopE



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.