



Tyrosine-protein phosphatase YopH (yopH)

OVERVIEW

PRODUCT REFERENCE	PX-P4379
ORIGIN SPECIES	
FRAGMENT	Met12-Ser399
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	52.10kDa
FORM	
BUFFER	20 mM Tris-HCl pH8.0, 100 mM NaCl, 1 mM DTT, 0,1% sodium azide, 0.1mM EDTA, 10% glycerol
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P15273
ALIAS/SYNONYMS	Tyrosine-protein phosphatase YopH, Virulence protein

