



Plastid acetyl-CoA carboxylase (Acc-1)

OVERVIEW

PRODUCT REFERENCE	PX-P4383
ORIGIN SPECIES	
FRAGMENT	Met1479-Lys2311
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	93.74kDa
FORM	
BUFFER	PBS pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	B2ZGL3
ALIAS/SYNONYMS	Plastid acetyl-CoA carboxylase

