



hSAss-TwinStrep-MYBPHL_iso (MYBPHL)

OVERVIEW

PRODUCT REFERENCE	PX-P4387
ORIGIN SPECIES	
FRAGMENT	Met1-Asn354
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	41.34KDa
FORM	
BUFFER	50mM Tris-Cl pH 7.5, 150mM NaCl
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	A2RUH7
ALIAS/SYNONYMS	MYBPHL/Myosin-binding protein H-like

