



Protein YopD(yopD)

OVERVIEW

PRODUCT REFERENCE	PX-P4403
ORIGIN SPECIES	
FRAGMENT	Ile61-Val306
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	28.01KDa
FORM	
BUFFER	TBS pH 7.5, 4 M urea
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P37132
ALIAS/SYNONYMS	yopD/Protein YopD

