



Mucin-4(MUC4) Nter GST

OVERVIEW

PRODUCT REFERENCE	PX-P4406
ORIGIN SPECIES	
FRAGMENT	Gln1875-Phe1914
TAG	N-terminal Gst Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	30.75kDa
FORM	
BUFFER	50 mM Tris-HCl , pH 8.0, 150 mM NaCl
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q99102
ALIAS/SYNONYMS	MUC4/Mucin-4

